

Product datasheet for **RR209530L2V**

Wisp2 (NM_031590) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wisp2 (NM_031590) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Wisp2
Synonyms:	CCN5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_031590
ORF Size:	750 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209530).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031590.1 , NP_113778.1
RefSeq Size:	1741 bp
RefSeq ORF:	753 bp
Locus ID:	29576
UniProt ID:	Q9JHC6
Cytogenetics:	3q42
Gene Summary:	involved in negative regulation of cell growth; mRNA expression is lost upon cell transformation induced by activated H-ras oncogene and inactivated p53 tumor supressor [RGD, Feb 2006]



[View online »](#)