

Product datasheet for **RR209551L4V**

Tas2r118 (NM_023994) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tas2r118 (NM_023994) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tas2r118
Synonyms:	T2R3; Tas2r16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023994
ORF Size:	897 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209551).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023994.1 , NP_076484.1
RefSeq Size:	900 bp
RefSeq ORF:	900 bp
Locus ID:	78980
UniProt ID:	Q9JKU0
Cytogenetics:	4q22
Gene Summary:	member of a family of mammalian taste receptors that couple with the G protein alpha subunit gustducin; may play a role in bitter taste perception [RGD, Feb 2006]


[View online »](#)