

## Product datasheet for **RR209557**

### Chst5 (NM\_001107431) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Chst5 (NM\_001107431) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Chst5  
**Synonyms:** RGD1561144  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209557 representing NM\_001107431  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGCCCCGGGCGGCAGCGCAGAGCAGCACAGAAGCGCAGCCGGCAGCGCAGGCTCTAGCCAGAGGCA  
TGGCGCTACCCCGCTTCCAGCACTGTAGTGCTTCGCTCCTGATAGTGCAGACCGGCATCCTGGTCTT  
CCTGGTCTCCCGCAAATGCCATCATCCCAGCAGGCCGTGGGAGCATGTGCATGTGCTGGTGCTGCC  
TCCTGGCGCTCGGCTCGTCTTTCGTGGCCAGCTCTTACGCCAACACCCGGATGTCTTCTACCTGATGG  
AGCCAGCTTGGCAGTGTGGGATGCGTTGTGCGAGGCCAGTGCCCCGCGCTCCACATGGCCGTGGCGCA  
CCTGATCCGTTCACTGTTCTATGCGACATGGACGTGTTGATGCTACCTGCCCTGGCGCCGAATATC  
TCGGATCTCTCCAGTGGCAGTGAGCCGGCGTTGTGCTCGCCTCCAGTCTGCGACGCCTTCGCTCGGG  
GCAACATTAGCAGCGAGGAGGTGTGAAGCCTCTGTGTGCAACCGGCCCTTCGGCTGGCTCAGGAAGC  
CTGCAGCTCTACAGTACGTAGTGCTCAAGGAGGTGCGCTTCTTAACTGCAGGTGCTCTACCCCTG  
CTCAGCGACCTGCGCTCAACTGCGCATCGTGCACCTAGTACGCGACCCGCGGGCCGTGCTACGTTCCC  
GAGAGCAGACAGCAAGGCGCTGGCAGGGACAATGGCATGTCCTGGGTACCAACGGCAGCATGGGTGGA  
GGCAGACCCCGGCTGCGGTGGTCAAGGAGGTGTGCCGAGCCATGTGCGCATCGCAGAGGTAGCATTG  
CACAAGCCGCGCCTTCTTGAAGATCGTACCCTGTTGCGCTATGAAGACCTGGCCCGAGACCCGC  
TCACCGTAATCCGTGAATCTATGCCTTACTGGTCTGGGTCTCAGCCGAGCTCCAGACTTGGATTTA  
CAACATCACCCATGGCTCGGGCCAGGCGCAGCCGTGAGGCCTTCAAGACCACATCCAGGGATGCGCTC  
AGTGTATCCAGGCTTGGCGCCACACACTGCCCTTGGCAAGATTGCGCGGTCCAGGAGCTGTGCGCGG  
GTGCACTGCAGTGTGCGTTACCGCCTGTGCACTCGGAGCTTGAACAACGGACCTCTCTGGAACCT  
CCTGCTGCCAGGCGATGGACAGTTTCAAGTGGCGTGTCCACGGATAAGCAACCGGAATCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR209557 representing NM\_001107431  
Red=Cloning site Green=Tags(s)

MRPGRAAQSSTEAPAAQALARGMRLPRFSSTVLSLLIVQTGILVFLVSRQMPSSPAGRGEHVHVLVLS  
 SWRSGSSFVQQLFSQHPDVFYLMPEAWHVWDALSQASAPALHMAVRDLIRSVFLCDMDVFDAYLPWRRNI  
 SDLFQWAVSRALCSPVCDAFARGNISSEEVCKPLCATRPFGLAQEACSSYSHVVLKEVRFNQLVLYPL  
 LSDPALNLRIVHLVRDPRAVLRREQTAKALARDNNGIVLGTNGTWVEADPRLRVSEVCRSHVRIAEVAL  
 HKPPPFLQDRYRLVRYEDLARDPLTVIRELYAFTGLGLTPQLQTIWYNIITHGSGPGARREAFKTTSRDAL  
 SVSQAWRHHTLPFAKIRRVQELCAGALQLLGYRPVHSELEQRDLSLDLLLPRGMDSFKWASSTDKQPES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001107431

<b>ORF Size:</b>	1254 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001107431.1</a> , <a href="#">NP_001100901.1</a>
<b>RefSeq Size:</b>	1991 bp
<b>RefSeq ORF:</b>	1257 bp
<b>Locus ID:</b>	307859
<b>Cytogenetics:</b>	19q12
<b>MW:</b>	46.9 kDa