

## Product datasheet for **RR209569**

### **Arl13b (NM\_001107101) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Arl13b (NM\_001107101) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Arl13b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209569 representing NM\_001107101  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGTTCAATCTGATGGCCAACCTGCTGCAACTTGTTC AAGCGGTGGCGGGAGCCTGT CAGAAAGGTGACAC  
 TTGTGATGGTTCGACTTGATAATGCTGGTAAAACAGCCACAGCAAAGGGAATCCAAGGAGAGCATCCTGA  
 AGACGTGGCTCCCACTGTTGGATTTTCTAAAATTGATCTGAGACAAGGAAAGTTTGAAGTTACCATCTTT  
 GACTTAGGAGGTGAAAAAGAATTCGGGAATTTGGAAGAATTATTATGCTGAATCCTATGGGGTAAATAT  
 TTGTTGTGGATTCCAGTGATGAAGACAGAATGGAAGAAACGAAGGAGACAATGTCAGAAGTGCTAAGACA  
 CCCTCGGATATCAGGAAAGCCTATATTGGTGTGGCAAATAACAAGACAAGGAGGGGGCTTTAGGAGAA  
 GCTGATGTGATTGAGTGTCTGTCTCGGAAAAGCTGGTCAATGAGCACAAGTGCCTGTGTCAGATAGAAC  
 CTTGTT CAGCAGTCTTGGGATATGGAAGAAAATCGACAAGTCCATTA AAAAGGGGCTTTATTGGCTACT  
 GCATATCATTGCAAAGATTTTGTGCTTAAGTGAACGCATCCAGAAAGACACAACCTGAACAGCGAGCT  
 CTTGAGGAACAAGAGAAATGTGAGAGGGCTGAAAGAGTAAGGAAGTTACGAGAAGAAAGAGAACAAGAAC  
 AGACGGAACCTTGATGGGACCAGTGGTCTGGCAGAGATGGACTCAGGACCAGTTCTTACAAATCCTTTCCA  
 ACCCATAGCAGCTGTATCATTGAGAATGAAAAGAAACAGGAAAAAGAGAAAAAGAAGACTGCGGAG  
 AAGGACAGTGATGTTGGCCACCCAGAACATAAAGGGGAACAGGAGGGAGCAGAGTCCGAGATGAGGCGAG  
 GTGGCTGCCTCAGGAACCTCACAGGGATGTGGTAAATAGTTACAAGGAGGCACTGT CACAGCAACTGGA  
 CAACGAGGATGAGCTAGACCAACGAGTGT CAGAACCAGGTGATAACAGTAAGAAGAAAACCAAGAAATTA  
 AGAATGAAAAGGAGTCATCGGGTAGAACAGTCAATACAGACGAGTCTGCTCCTAAGAGTCCCACGCCGC  
 CCCACCTCCCCCTCCTGTTGGCTGGGGAACCCCAAAGTCACTAGACTTCCAAAACCTTGAGCCTCTTGG  
 TGAACACGTCATAATGATTTCTATGGGAAGCCGCTGCCTCCCTGGCTGTGCGACAAAGACCTAATGGT  
 GATGCTCAGGACATGATCTCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



**Protein Sequence:** >RR209569 representing NM\_001107101  
Red=Cloning site Green=Tags(s)


MFNLMANCCNLFKRWREPVRKVTLMVMGLDNAGKTATAKGIQGEHPEDVAPTGVGFSKIDLRQGKFEVTIF  
 DLGGGKRIRGIWKNYYAESYGVIFVVDSSDEDRMEETKETMSEVLRHPRI SGKPILVLANKQDKEGALGE  
 ADVIECL SLEKL VNEHKCLCQIEPCSAVLGYGKKIDKSIKGL YWLLHIIAKDFDALSERIQKDTTEQRA  
 LEEQEKCERAERVRKLEEREQEQTLDGTSGLAEMDSGPVLTNPFQPIAAVI IENEKKQEKEKKKKTAE  
 KSDSDVGHPEHKGEQEGAESQNEAGGCLRNPHRDVMNSYKEALSQQLDNEDEL DQRVSEPGDNSKKKTKKL  
 RMKRSHRVEPVNTDESAPKSPTPPPPPVWGTPKVTRL PKLEPLGETRHNDFYGKPLPPLAVRQRPNG  
 DAQDMIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001107101
<b>ORF Size:</b>	1281 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001107101.1</a> , <a href="#">NP_001100571.1</a>
<b>RefSeq Size:</b>	3535 bp
<b>RefSeq ORF:</b>	1284 bp
<b>Locus ID:</b>	304037
<b>Cytogenetics:</b>	7q11
<b>MW:</b>	48 kDa