

Product datasheet for **RR209586**

Prr22 (NM_001126375) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prr22 (NM_001126375) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Prr22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR209586 representing NM_001126375 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACTGGTAAACAGCCAGTTTTCAGGGCTCCTCAGAGGCCAAGGCTGTGTCTATGTCCTGGGCTTTGG
GACCTCAACCCCATGGCCTGGTTCCCCATGGAGCCCCACTCAGGTGTCCACCACCCACAGGCTTCCA
GATGGCTCCCTGCGGCTGCTTCTTCGACCCCGCATCTACCGCATCGAGTGGGCCACCTCCGACTTTGGC
CAGTCATCCCTGTACAAGGTGGCAGTGGCTGGGGTCTGCCTTATCTGGTGGCTACCTGCTGGAGGCCC
CGTCCTACCTCAAGACCCCGGGGCCACCACTCTGCTTTATCCCCACTACCAGCCGGCCCCAGCGGGCC
ACAGTACCTCACACATTTCTTCCACCAGAGGAGCCTATGCCTGAGGCTCTGGGCTTCATGCGGGACGGA
GGGCCCCTCAACTTCATGGAGGTGCTCAGAGACGGCTGGCACCCCGCCAACCCCAAAGAGACTAAGC
CATCTCCCTTGCTTATCAGAGTCCCCACAGCACACAGCTGCCACCTGGGCCCTATGGCCATCTCGGTGG
CCAGACCAGCCAGTTTCCAGGACCCAGGTGACCATGAGGCCACCGAGGCCCCAGGGAGCTGCAGAGT
AGTGGTGCAGCCAGGCCGGGTCTGCAGTTCCCTCCTGGGCTGGAGAGCCACGGTGGCAAGGTGGAGG
AGGTCACCCAGTGGGCTCAAGTGAGACCGTGCCTCCGAGGTGGCTCGTGCCTTCTTCTACCAGACAA
GGTCCTTCTAGAAGATGCCATGAACTTTTGTGTTCTCCTGGTGGCACTGAACCTGAGGTGGCCTTG
CACAGGGTCTTGCCCCGAGTGGGGACAGTGGTGGTGGGGACGACTTCCCAACTGACATCCGCT
CCCTGCACCTGCCAGATGACCTGCTGTCTTTGACTACAGTGTCCCTGAGATTCTGGACGCGGTGGCCAA
CGTGGATTATGTTTTAGCTTCAAGGCCCTGGAGGATGAGCCGCTGCCCATCTTGGGGCCCTGTCACT
GACACGGCGGCCCTGGCCTGCGGTCCCATCAACCGGGGAAAAACCGAGCACATCCACCAAGAAAGGA
AACCAAGCGGTAGACACAGGCAGGCCACTGGCCTTGCTGGCACCCTGCCACAGGACCCAGGCTGGACCC
AGGAGCCACCCCAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »


Protein Sequence: >RR209586 representing NM_001126375
 Red=Cloning site Green=Tags(s)

MTGKQPVFKAPQRPRLCPMSWALGPSTPMWFPHGAPLRCPPPTGFQMAPCGCFFDPRIYRIEWATSDFG
 QSSLYKVAVAGGPALSGGYLLEAPSYLKTGPPPLLYPHYQPAPSGPQYLTHFLPPEEPMPALGFMRDG
 GPLNFMVLRDGLAPPTPKETKPSLLITVPTAHTLPPGPYHGLGGQTSQFPGPQVTMRPTEAPRELQS
 SGAARPLQFPPGPGEPTVAKVEEVTVPGSSETVPPEVARAFFLPDKVLLEDAMKLFDCPLPGGTEPEVAL
 HRGPGPLRDSGGGGDDFPTDIRSLHLPDDLSSFDSVPEILDVANVDYVFSFKALEDEPLPHLGAPVT
 DTAAPGLRSHQPGKKPSTSTKKGKPSGRHRQATGLAGTTATGPRLDPGATPN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001126375

ORF Size: 1206 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001126375.1, NP_001119847.2

RefSeq Size: 1766 bp


RefSeq ORF: 1209 bp

Locus ID: 680177

Cytogenetics: 9q12

MW: 42.6 kDa

Product images:



Circular map for RR209586