

Product datasheet for **RR209595**

Ttc30b (NM_001127607) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ttc30b (NM_001127607) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ttc30b
Synonyms:	MGC188467
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR209595 representing NM_001127607
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCAGGGCTGAGCAGCTCTCAGATCCCCGACGGTGAGTTCACAGCGGTGGTGTACAGGCTCATCCGCG
ACTCCCCTACTCAGAGCGGTGACGCTGCTGAGCGCAGAGCTGCAGGGAGCCCCGCAGCCGCGCAGG
GCTGTCGCTGCTGGCCTACTGCTACTACCGCTTGCAGGAGTTCGAGCTCGCTGCGGAGTGTATGAGCAA
CTGAGCCAGATGCATCCGGAACCGAGCAGTACCGCTGTACCAGGCCAGGCGCTCTACAAGGCCTGCC
TGTATCCAGAGGCCACCCGGTTCGCTTCTCTGACAAACCCACCTATCAGACTCGAGTCTTCTGCTCT
GCAGGCTGCTATCAAGTACAGTGAAGGGGATCTCCAGGAGCCAGGAGCCTAGTGGAGCAGCTGCTGAGT
GGGAAGGAGCTGAGGACAGTGGAGGGGAGAATGATTATGATGGTCAGATCAACCTGGTGTCTGCTCT
ACAAGGAGGGACACTATGAAGCTGCCTGTTCCAAGTCTTAGCGGCCCTGCAGGCTTCTGGTTACCAGCC
TGACATTTCTACAACCTGGCTTTAGCCTATTACAGCAGTCGGCAGTATGCGCCTGCTCTGAAGCATATC
GCTGACATCATTGAGCGTGGCATCCGTCAACCCAGAAGTGGTGGGATGACCACCGAAGGCATCG
ATGTTTCAAGTGTGGCAACACCATAGTCTTCCAGACTGCTCTGGTGAAGCCTTCAACCTCAAGGC
AGCCATAGAGTACCAGCTGAGAACTATGAAGCAGCTCAGGAAGCCCTCACTGACATGCCACCCAGAGCA
GAGGAAGAGTTGGACCCCGTACCCCTGCACAACCCAGGCTCTGATGAACATGGATGCCAAGCCACGGAGG
GCTTTGAAAAGCTCCAGTTTCTGCTCCAGCAAAACCCCTCCCCCAGAGACCTTTGGCAACCTGCTGCT
GCTTTACTGTAATATGAGTATTTGACCTGGCAGCTGATGTTCTAGCAGAAAACGCCACTTGACTTAC
AAGTTCCTCACACCCTATCTCTATGACTTCTGGATGCCATGATCACTTGCCAGACAGCTCCTGAAGAGG
CTTTCATTAACCTTGATGGGCTAGCAGGCATGCTGACGGAGCAGCTCAGGAGACTCACCATAAAGTACA
GGATGCAAGACATAGTAGAGATGATGAATCTGCCAAAAAGGCTGTAACGATTATGATGAAACCTTGAG
AAGTATATCCCTGTACTGATGGCCAGGCCAAAATCTACTGGAATCTTGAGAATTATCCAATGGTGGAGA
AGATCTTCCGCAATCTGTGGAATTCTGTAATGACCACGATGTGTGGAAGCTGAATGTGGCCATGTCTCT
CTTCATGCAGGAAAACAAGTACAAGGAGGCCATTGGTTTCTATGAGCCATTGTCAAGAAGAATTATGAT
AACATCTGAGTGTGACGCTATCGTATTAGCTAACCTCTGTGTTTCTACATTATGACAAGTCAAAATG
AGGAAGCTGAGGAGCTGATGAGGAAGATTGAGAAGGAGGAGGAACAACCTCCTATGATGACCCGGACAA
GAAAATCTACCACCTTTGCATTGTGAATTTGGTATAGGAACCTTTATGTGCCAAAGGAACTATGAC
TTTGGCATCTCTCGGTTATCAAGAGCCTGGAGCCTTACCATAAGAAGCTGGGAAGTATGACGTGGTATT
ATGCCAAAAGATGCTTCTGTCCTTACTGGAGAACATGTCAAAGCACACAATCATGCTGAGGACAGCGT
CATTCAAGAATGTGTCCAGTTTCTAGAGCACTGTGAACCTTATGGCAGAAGCATCCCTGCTATTATCGAA
CAGCCCTTGGAGGAGGAAAGAATGCACACGGGGAAGAATACGGTCACTTATGAGTCCAGAAAATTAAGGG
CTCTGATTTATGAGATCATTGGGTGGAATATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR209595 representing NM_001127607
 Red=Cloning site Green=Tags(s)


MAGLSSSQIPDGEFTAVVYRLIRDSRYSEAVQLLSAELQGSPPSRAGLSLLAYCYRRLQEFELAAECYEQ
 LSQMHPELEQYRLYQAQALYKACLYPEATRVAFLLDNPTYQTRVLRQAQAIKYSEGDLPGARSLVEQLLS
 GEGAEDSGGENDYDQGQINLGLLLYKEGHYEAACSKFLAALQASGYQPDISYNLALAYSSRQYAPALKHI
 ADIIERGIRQHPQLGVGMTTEGIDVRSVGNITIVLHQTALVEAFNLKAAIEYQLRNYEAAQEALTDMPRA
 EEELDPVTLHNQALMNMDAKPTEGFELQFLQNPFPPEFTGNLLLLYCKYEYFLAADVLAENAHLY
 KFLTPYLYDFLDAMITCQTAPEEAFIKLDGLAGMLTEQLRRLTIQVQDARHSRDESAKKAVNDYDETLE
 KYIPVLMQAQAIYWNLENYPMVEKIFRKSVEFCNDHDVWKLNVAVHVLFMQENKYKEAIGFYEPVVKKNYD
 NILSVSAIVLANLVCVSYIMTSQNEEAELMRKIEKEEQLSYDDPKKIYHLICIVNLVIGTLYCAKGNYP
 FGISRVIKSLEPYHKKLGTDTWYAKRCFLSLEENMSKHTIMLRDSVIQECVQFLEHCEL YGRSIPAIIE
 QPLEEERMHTGKNTVTYESRKL RALIYE IIGWNM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001127607

ORF Size: 1992 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127607.1](#), [NP_001121079.1](#)

RefSeq Size: 2589 bp

RefSeq ORF: 1995 bp

Locus ID: 499814

UniProt ID: [B2RYD6](#)

Cytogenetics: 3q23

MW: 76.1 kDa

Gene Summary: Required for polyglutamylation of axonemal tubulin. Plays a role in anterograde intraflagellar transport (IFT), the process by which cilia precursors are transported from the base of the cilium to the site of their incorporation at the tip.[UniProtKB/Swiss-Prot Function]