

## Product datasheet for **RR209628**

### Slc6a20 (NM\_133296) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Slc6a20 (NM_133296) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc6a20
Synonyms:	rB21a; Xtrp3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR209628 representing NM\_133296  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAGATTAGCAATTAAGGCGGGCGAGCCGCGCCAGAGACCAGGCCCTGACGAGAAGCGAGCGCGGG  
 ACATGGAGAAGGCACGGCCTCAATGGGGCAATCCGCTGCAGTTTCGTTTTGCGCTGTATCTCTACGCGT  
 GGGTTTGGCAATGTGTGGCGCTTCCCCTACCTGTGCCAGATGTACGGCGGAGGGAGTTTCTGGTCCCC  
 TACCTCATCATGCTCATTGTGGAGGGGATGCCACTCTTGTACCTGGAGCTGGCTGTGGGGCAGCGCATGC  
 GGCAGGGCAGCATTGGTGCCTGGAGGACCATCAGCCCCTACCTAGTGGTGTGGTGTGCCAGTGGT  
 GGTCTCCTTCTCCTCTCCATGTACTACAATGTCATCAATGCCTGGGGCTTCTGGTACCTCTTCCACTCC  
 TTCCAGGATCCCCTGCCGTGGTCTGTCTGCCACTGAATAGTAACCGCACAGGCTATGATGAGGAGTGTG  
 AGAAGGCTTCGTCGACACAGTACTTCTGGTACAGGAAAACACTCAACATCTCACCGTCCATCCAGGAGAA  
 TGGAGGAGTGCAGTGGGAGCCAGCCCTGTCCCTACCCTGGCCTGGCTGATGGTATATCTGTGCATCCTG  
 AGAGGCACCGAATCTACTGGCAAGGTGGTCTACTTCACCGCATTGATGCCTACTGTGTTCTTATTATCT  
 ACTTGGTCCGTGGCCTCACACTCCATGGAGCCACCAATGGCCTGATGTACATGTTACACCTAAGATTGA  
 GCAGCTAGCCAACCCCAAGGCTGGATCAATGCAGCCACGCAGATCTTCTTCTACTGGGCTTGGGTTTT  
 GGCAGCCTGATCGCTTTTGGCAGCTACAATGAACCCTCCAACGACTGCCAGAAGCATGCTGTATTGTGT  
 CTGTCAACAACAGCTCCACCTCCATATTTGCCAGCATTGTACCTTCTCCATCTATGGCTCAAGGCCAC  
 CTTCAACTATGAAAACGCTTAAACAAGGTGATTCTGTGCTGACCAATCTTTTGACCTGAAGATGGC  
 TTTCTGACAGCCAGCAACCTGGAGGAGGTGAAGGACTACCTGGCATCTACTTACCCAAACAAGTACAGTG  
 AAGTGTTCACACATTAGAACTGCAGCTTGAATCAGAGCTGAACACGGCTGTCCAAGGCACAGGCCT  
 GGCTTTCATCGTCTACGCTGAGGCCATTAACAACATGGAAGTGTCCAGCTCTGGTCACTGCTCTACTTC  
 TTCATGCTGCTGATGCTGGGAATGGGGAGCATGCTTGAATACAGCGGCCATCCTCACCCCTCTGACTG  
 ACAGCAAGGTCATCTCCAGCTACCTGCCAAGGAGGCCATTTCAAGTCTGGTGTGCCTCATTAACTGTGC  
 TGTTGGCATGGTGTCCACATGGAGGCTGGGAACACTGGTTTGACATATTCAATGACTATGCAGCCACG  
 CTGTCTCTGCTGCTCATTGTGCTGGTGGAGACTATAGCTGTGTGCTACGTGTATGGGCTGAGGAGATTTG  
 AAAGTGTCTTCGGGCCATGACTGGCCGGCCCTCACTGGTACTGGAAGGCCATGTGGGCTTTTGTGAG  
 CCCACTGCTCATCATCGGCTCTTATCTTCTACCTGAGTACTACATCCTCACGGGAACGCTGCAGTAC  
 CAAGCCTGGGATGCTACTCAGGGCAGCTGGTGACCAAGGATTACCTCCACATGCACTAGCTGTATCG  
 GTTTGCTGGTGGCTTCTACTATGTGCATCCCCTGGTGGCCTGGGGACTTTCATCAGGAATCGCCT  
 CAAGAGGGGAGGCTTCTCCCACTGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>RR209628 representing NM\_133296  
 Red=Cloning site Green=Tags(s)

MRLAIKRASRGQRPGPDEKRRARMEKARPQWGNPLQFVFIACISYAVGLGNVWRFPPYLQMYGGGSLVLP  
 YLIMLIVEGMPLLYLELAVGQRMQRQGSIGAWRTISPYLSGVGVASVVVFFLSMYYNVINAWGFWYLFHS  
 FQDPLPWSVCPNLSNRTGYDEECEKASSTQYFWYRKTNLISPSIQENGGVQWEPALCLTLAWLMVYLCIL  
 RGTESTGKVVYFTALMPYCVLIIYLVRLTLHGATNGLMYMFTPKIEQLANPKAWINAATQIFFSLGLGF  
 GSLIAFASYNPSNDCQKHAVIVSVINSSTISFASIVTFSIYGFKATFNENCLNKVILLTNSFDLEDG  
 FLTASNLEEVKDYLASTYPNKYSEVFPHIRNCSLESELNTAVQGTGLAFIVYAEAIKNMEVSQLWSVLYF  
 FMLLMLMGMSMLGNTAAILTPLTDSKVISSYLPKEAISGLVCLINCAVGMVFTMEAGNYWFDIFNDYAAT  
 LSLLLIVLVETIAVCYVYGLRRFESDLRAMTGRPLNWWKAMWAFVSPLLIIIGLFIYFVSDYILTGLTQY  
 QAWDATQGLVTKDYPPHALAVIGLLVASSTMCIPLVALGTFIRNRLKRGSSPVA


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_133296

ORF Size: 1848 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_133296.1</a> , <a href="#">NP_579830.1</a>
<b>RefSeq Size:</b>	1950 bp
<b>RefSeq ORF:</b>	1851 bp
<b>Locus ID:</b>	113918
<b>UniProt ID:</b>	<a href="#">Q64093</a>
<b>Cytogenetics:</b>	8q32
<b>MW:</b>	69 kDa
<b>Gene Summary:</b>	may act as a neurotransmitter transporter [RGD, Feb 2006]