

## Product datasheet for **RR209630**

### **Fbxo22 (NM\_001037770) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fbxo22 (NM\_001037770) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Fbxo22  
**Synonyms:** MGC124731  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209630 representing NM\_001037770  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAGCCAGCAGGCGGGCTCCTCCTCCACAGATCCGCGAGGCACCTATGTGCTGAGTAACCTTGCGG  
AGGTGGTGGAGCGTGTGTTACCTTCTGCCGGCAAAGCGTTGCTTCGGGTAGCCGGAGTATGCCGCCT  
GTGGAGAGAGTGTGTGCGTAGAGTGCTGCGGACCCATCGCAGCGTGACCTGGATCTCCGCGGGTGTGGCG  
GAGGCGGGCCACCTGGAGGGACATTGCTTGGTCCGCGTGGTAGCTGAAGCGCTTGAGAATGTCCGCATCT  
TACCACGGACAGTTCTCTACATGGCAGATTCTGAAACGTTTCATCAGCCTGGAGGAGTGCCGGGGCCATAA  
AAGAGCAAGGAAAAGAACTACTATGGAGACAGCATGTGCCTGGAGAAGCTGTTCCCAAACAGTGTGAG  
GTCCTTGGGATTGTGACCCAGGAATTGTAGTGACTCCAATGGGATCAGGTAGCAATCGGCCTCAGGAAA  
TAGAAATTGGAGAATCTGGTTTTGCTTTATTATCCCTCAAATGAAGGAATAAAAATTC AACCCTTTCA  
TTTTATTAAGGACCCCAAGAATTAACCTCTCGAAAGACATCAACTACTGAAGTGGGTCTTCTGGATAAC  
CCCGAACTCCGTGTGGTCTTGTCTTTGGCTATAACTGCTGT AAGGTGGGAGCCAGTAATTACCTGCACC  
GGGTGGTCAGCACTTTCAGCAACATGAACATCATCTTGGCTGGAGGCCAGTGGACAATGTCTTCACT  
GACTTGTGAGAAGAACCCTCTGGATATTGACGCCACAGGTGTGGTCCGACTGTCATTTAGTGGCACCGA  
ATCCAGAGTGCCACGGTCTCCTCAACGAGGATGTAATGATACCAAGACCGCTGAGGCCGCCATGCAGC  
GCCTCAAAGCAGCCAACATCCCTGAACAGAACACCATTGGCTTCATGTTTGCCTGTGTTGGCCGAGGCTT  
TCAGTACTACAGAGCCAAGGGGAATGTTGAGGCTGATGCATTTAGAAAAGTTTTTCCCAATGTCCCTTA  
TTTGGCTCTTTGGAAATGGGGAGATCGGCTGCGACCGGATAGTCACTGGGAACTTTATCCTGAGGAGGT  
GTAATGAGGTGAAGGAGGAGGATCTGTTCCATAGCTACACAACCATCATGGCTCTCGTTCACCTGGGGC  
CTCCAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR209630 representing NM\_001037770  
Red=Cloning site Green=Tags(s)

MEPAGGSSSTDPRTGYVLSNLAEVVERVFTFLPAKALLRVAGVCRLWRECVRRVLRTHRSVTWISAGVA  
 EAGHLEGHCLVRVVAEALENVRLPRTVLYMADSETFISLEECRGGHKRARKRTTETACALEKLFKQCQ  
 VLGIVTPGIVVTPMGSGSNRPQEIIEIGESGFALLFPQIEGIKIQPFHFIKDPKNLTLEHQLTEVGLLDN  
 PELRVVLFVGYNCCKVKGASNYLHRVVSTFSNMNIILAGGQVDNLSSLTCEKNPLDIDATGVVGLSFSGHR  
 IQSATVLLNEDVNDTKTAEAAMQRLKAANIPEQNTIGFMFACVGRGFQYYRAKGNVEADAFRKFPPNVPL  
 FGFFGNGEIGCDRIVTGNFILRRCNEVKEEDLFHSYTTIMALVHLGASK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001037770

<b>ORF Size:</b>	1197 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001037770.1</a> , <a href="#">NP_001032859.1</a>
<b>RefSeq Size:</b>	1832 bp
<b>RefSeq ORF:</b>	1200 bp
<b>Locus ID:</b>	300724
<b>Cytogenetics:</b>	8q24
<b>MW:</b>	44 kDa