

## Product datasheet for **RR209633L4V**

### Dll1 (NM\_032063) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dll1 (NM_032063) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dll1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032063
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032063.2</a> , <a href="#">NP_114452.1</a>
RefSeq Size:	3171 bp
RefSeq ORF:	2145 bp
Locus ID:	84010
UniProt ID:	<a href="#">P97677</a>
Cytogenetics:	1q12
Gene Summary:	homolog of the Notch Delta ligand; member of the delta/serrate/jagged family; may play a role in cell fate decisions during hematopoiesis [RGD, Feb 2006]



[View online »](#)