

Product datasheet for **RR209671L4V**

Hnrnph2 (NM_001014019) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnph2 (NM_001014019) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnph2
Synonyms:	Hnrph2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001014019
ORF Size:	1347 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209671).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001014019.1 , NP_001014041.1
RefSeq Size:	1725 bp
RefSeq ORF:	1350 bp
Locus ID:	308650
UniProt ID:	Q6AY09
Cytogenetics:	Xq32



[View online »](#)

Gene Summary:

This protein is a component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes which provide the substrate for the processing events that pre-mRNAs undergo before becoming functional, translatable mRNAs in the cytoplasm. Binds poly(RG) (By similarity).[UniProtKB/Swiss-Prot Function]