

## Product datasheet for **RR209673L3V**

### **Pdgfc (NM\_031317) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	Pdgfc (NM_031317) Rat Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	Pdgfc
<b>Synonyms:</b>	SCDGF
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
<b>Tag:</b>	Myc-DDK
<b>ACCN:</b>	NM_031317
<b>ORF Size:</b>	1035 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RR209673).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_031317.1</a> , <a href="#">NP_112607.1</a>
<b>RefSeq Size:</b>	1116 bp
<b>RefSeq ORF:</b>	1038 bp
<b>Locus ID:</b>	79429
<b>UniProt ID:</b>	<a href="#">Q9EQX6</a>
<b>Cytogenetics:</b>	2q33-q34
<b>Gene Summary:</b>	expressed in the embryonic spinal cord floor plate and ventricular zone, adult expression is localized to motoneurons in several regions of the brain [RGD, Feb 2006]



[View online »](#)