

Product datasheet for **RR209698**

Spata3 (NM_001108805) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spata3 (NM_001108805) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Spata3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209698 representing NM_001108805
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAAGAAGGTCAAAAAGAAAAAGTCTGATTCCAGGCGCCGAGGAATCCATTTACCCGAGACTAGCT
 CCGACTGCTCACACAGGCCAGCTCAGAAACCCAGCGGTATGTCAGACCCCTGTTACCACCCGAGAC
 ACAGTCTCTCCAGGAGTCCATCCACCACCAAATTAGCCCTGAGCCTAAACTGCAACAGCCTTCATCA
 CAACCTCTCTACCTTCAGAAGAAGTGTCTCATCTCCGTGCCTGGTATCCACGGAACCTAAGCCTCCAT
 GCAGAAAAGCAGCAGTGCCTCTGACCTATGCGGCCCCACGCTCCTGCTCCTGTGTGGTCTGCCCTGGCAG
 CTCCGCTTGTGGCGTCTGTTGGGCTTTGCCATAGCCGAATCTTTGATGTCTTCTGCCTCGGGACTGG
 TCATCCATGCCAGGGAGAGGATTCCCAAACCTTCTCACCTTCTACAGAAAACCTTCAAGGAAATACTGCA
 CTTCCCGTAACTCTCGTCTTCAAGCCCTCGGAAGTGTGCTGTGGCTCTGGGGCCCTGGGAGTTGCT
 ACTTCATGGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR209698 representing NM_001108805
 Red=Cloning site Green=Tags(s)

MKKVKKKSDSRRRRNSISPQTSSDCSHRPSSETQRSCPDPSPPPQTQSLQESIPPHQISPEPKLQPPSS
 QPLLPSSEERVSSPCLVSTEPKPPCRKAAVPLTYAAPRSCSCVCPGSSACWRRRLGLCHSRIFDVLLPRDW
 SSMPPGRGFPNLLTFYRKPSRKYCTSRNSRASSPRNCCCGSGGPGSCLLHG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108805

ORF Size: 570 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108805.1, NP_001102275.1</u>
RefSeq Size:	1295 bp
RefSeq ORF:	573 bp
Locus ID:	363270
Cytogenetics:	9q35
MW:	20.9 kDa