

Product datasheet for **RR209706**

Elac1 (NM_001107406) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Elac1 (NM_001107406) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Elac1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR209706 representing NM_001107406 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTATGGATGTGACCTTTCTTGGGACAGGAGCAGCATATCCATCGCCCACTCGGGGCGCCTCTGCCG
TGGTTCTTCGGTGTGAGGGCGAGTGTGGCTCTTCGACTGTGGGAGGGCACACAGACCCAGCTTATGAA
AAGCCAGCTTAAAGCAGGGAGAATCACCAAGATCTTCATCACACATCTCCATGGAGACCACTTCTTTGGC
CTTCCAGGCCTCCTCTGCACCATCAGCCTGCAGAGTGGCTCTGTGGTCAGCAGGCAGCCTATTGAGATCT
ACGGCCCCGTGGCCTTCGGGACTTTATCTGGAGAACCATGGAGCTCTCTCACACGGAGTTGGTTTTCCC
TTATGTGGTCCACGAGCTGGTTCCCTACAGCTGACCAAGTGCCTGTGGAAGAACTAAGAGAATCTGCACAA
ATGGATGAAGCAGACAGCTCTCCCAAAGGACAAGGACGGAATCTTGCTAGATGCAGAGGAAAACCTCT
ACTGTCTGGTGGATGACGAACAGTTTGTGGTAAAAGCCTTCGCTCTTCCATAGAATTCATCTTTTGG
CTTTTCGCTCGTGGAAAAGAAGCGCACCGGAAAACCTCAGTGCTCAGAAGCTGAGAGACCTTGGTGTCCCA
CCAGGTCTGCCTATGGGAAGCTGAAGAATGGGATTTCTGTTGTTTGGATAATGGAGTTACAATTTCTC
CCCAAGATGTCTTAAAAAACCTATGGTTGGAAGAAAAGTCTGCATCCTGGGCGACTGTTCCGGGTCGT
GGGTGATGGGGGAGTGAAGCTGTGCTTTGAGGCAGACCTGCTGATCCATGAGGCGACTCTGGATGACTCC
CAGATGGACAAGGCGAAGGAGCATGGCCACAGCACACCACAGATGGCAGCAACATTTGCAAAGCTGTGCC
GTGCAAAGCGCCTCGTCTGACTCACTTCAGCCAGAGGTACAAACCAACTGCCTTAGCCAGAGAAGGGGA
GGCAGACGGCATCACGGAAGTGAAGCAAGCCGAGTCTGTGCTGGAGCTCCAAGAAGTACTCTAGCG
GAAGATTTTATGGTGATTGGCATTCCGATCAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR209706 representing NM_001107406
 Red=Cloning site Green=Tags(s)

MSMDVTLGTGAAYPSPTRGASAVVLRCEGECWLFDCGEGTQTQLMKSQKAGRITKIFITHLHGDHFFG
 LPGLLCTISLQSGSVVSRQPIEYGPVGLRDFIWRMELSHTEL VFPYVVHEL VPTADQCPVEELRESAQ
 MDEADSSPKGQGRTILLDAEENSYCLVDDEQFVVKAFRLFHRIPSFSGFSVVEKKRTGKLSAQKLRDLGVP
 PGPAYGKLKNGISVVLNDNGVTISPQDVLKPMVGRKVCILGDCSGVVGDGGVKLCFEADLLIHEATLDDS
 QMDKAKEHGHSTPQMAATFAKLCRAKRLVLTHFSQRYKPTALAREGEADGITELKKQAESVLELQEVTLA
 EDFMVGIIPIKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001107406

ORF Size: 1086 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001107406.1, NP_001100876.1

RefSeq Size: 2836 bp


RefSeq ORF: 1089 bp

Locus ID: 307604

Cytogenetics: 18q12.2

MW: 39.7 kDa

Product images:



Circular map for RR209706