

Product datasheet for **RR209711L3V**

Pigm (NM_024144) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|----------------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Pigm (NM_024144) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Pigm |
| Synonyms: | GPI-MT-I |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_024144 |
| ORF Size: | 1269 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR209711). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_024144.1 , NP_077058.1 |
| RefSeq Size: | 3700 bp |
| RefSeq ORF: | 1272 bp |
| Locus ID: | 79112 |
| UniProt ID: | Q9EQY6 |
| Cytogenetics: | 13q24 |
| Gene Summary: | mannosyltransferase that catalyzes the transfer of the first mannose to glycosylphosphatidylinositol in GPI anchor biosynthesis; activity may be localized to the luminal side of the endoplasmic reticulum membrane [RGD, Feb 2006] |



[View online »](#)