

Product datasheet for **RR209712L4V**

Lck (NM_001100709) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lck (NM_001100709) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lck
Synonyms:	Lck1; Lcktkr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001100709
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209712).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001100709.1 , NP_001094179.1
RefSeq Size:	1978 bp
RefSeq ORF:	1530 bp
Locus ID:	313050
UniProt ID:	Q01621
Cytogenetics:	5q36
Gene Summary:	mouse homolog is a kinase involved in activation-induced T cell death [RGD, Feb 2006]



[View online »](#)