

Product datasheet for **RR209725L3V**

Lgmn (NM_022226) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lgmn (NM_022226) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lgmn
Synonyms:	Prsc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022226
ORF Size:	1305 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209725).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022226.2 , NP_071562.2
RefSeq Size:	1843 bp
RefSeq ORF:	1308 bp
Locus ID:	63865
UniProt ID:	Q9R0J8
Cytogenetics:	6q32
Gene Summary:	putative cysteine protease [RGD, Feb 2006]



[View online »](#)