

Product datasheet for **RR209735L3V**

Crhbp (NM_139183) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crhbp (NM_139183) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Crhbp
Synonyms:	Crfbp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_139183
ORF Size:	966 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209735).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139183.2 , NP_631922.1
RefSeq Size:	1636 bp
RefSeq ORF:	969 bp
Locus ID:	29625
Cytogenetics:	2q12
Gene Summary:	binds corticotropin releasing hormone and modulates CRH action; may act as a reservoir of endogenous CRH in the brain [RGD, Feb 2006]



[View online »](#)