

## Product datasheet for **RR209746L4V**

### Rack1 (NM\_130734) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rack1 (NM_130734) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rack1
Synonyms:	Gnb2l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130734
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209746).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_130734.2</a> , <a href="#">NP_570090.2</a>
RefSeq Size:	1209 bp
RefSeq ORF:	954 bp
Locus ID:	83427
UniProt ID:	<a href="#">P63245</a>
Cytogenetics:	10q21
Gene Summary:	may act as a receptor for activated protein kinase C [RGD, Feb 2006]



[View online »](#)