

Product datasheet for **RR209773L3V**

Has2 (NM_013153) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Has2 (NM_013153) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Has2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013153
ORF Size:	1656 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209773).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013153.1 , NP_037285.1
RefSeq Size:	4172 bp
RefSeq ORF:	1659 bp
Locus ID:	25694
UniProt ID:	O35776
Cytogenetics:	7q33
Gene Summary:	membrane protein that synthesizes hyaluronan, a high molecular weight polysaccharide that is a component of the extracellular matrix; may be important for epidermal proliferation, thickness and differentiation [RGD, Feb 2006]



[View online »](#)