

## Product datasheet for **RR209785L4V**

### PPP1R9A (NM\_053473) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppp1r9a (NM_053473) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ppp1r9a
Synonyms:	Neb1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053473
ORF Size:	3285 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209785).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053473.2</a> , <a href="#">NP_445925.1</a>
RefSeq Size:	3846 bp
RefSeq ORF:	3288 bp
Locus ID:	84685
UniProt ID:	<a href="#">O35867</a>
Cytogenetics:	4q13
Gene Summary:	binds to protein phosphatase type I and to actin filaments; involved in regulation of synaptic transmission [RGD, Feb 2006]



[View online »](#)