

## Product datasheet for **RR209815**

### Socs2 (NM\_058208) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Socs2 (NM\_058208) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Socs2  
**Synonyms:** Cish2; Socs-2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR209815 representing NM\_058208  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGACCCTGCGGTGCCTGGAGCCCTCTGGGAATGGCGGGACAGGACACGGAGCCAGTGGGGACAGCGG  
 GGTGCGCCGGAGGATCAGTCCCCGGAGGCGCGCTCTGGCGAAGGCCCTGCGTGAGCTCAGTCAAACAGG  
 ATGGTACTGGGAAGTATGACTGTTAATGAAGCCAAAGAGAAATTAAGAGGCGCCAGAAGGAACCTTC  
 TTGATTAGAGATAGCTCCCACTCAGACTACCTATTAACCATATCTGTTAAGACATCAGCCGGCCGACTA  
 ACCTGCGGATTGAGTACCAAGACGGAAATTCAGATTGGATTCTATCATATGTGTCAAGTCCAAGCTTAA  
 GCAGTTTGACAGCGTGGTTCATCTGATTGACTACTACGTTGAGATGTGCAAGGACAAACGGACAGGCCCT  
 GAAGCCCCCGGAATGGGACTGTTACCTTTACCTGACCAAACCTCTGTATACATCAGCACCGACGCTGC  
 AGCACTTCTGCCGACTCAGCATTAAAGTGTACCGGTACCATCCGGGGACTACCTTTACCAACAAGACT  
 GAAAGATTACTTGAAGAGTATAAATTCAGGTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR209815 representing NM\_058208  
 Red=Cloning site Green=Tags(s)

MTLRCLPSGNGADRTRSQWGTAGSPEDQSPEARLAKALRELSQTGWYWSMTVNEAKEKLKEAPEGTF  
 LIRDSSHSYLLTISVKTSAGPTNLRIEYQDGKFRLDSIICVKSKLKQFDSVVHLIDYYVQMCKDKRTGP  
 EAPRNGTVHLYLTKPLYTSAPTLQHFRLSINKCTGTIRGLPLPTRLKDYLEEYKFQV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_058208

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_058208.1](#), [NP\\_478115.1](#)

**RefSeq Size:** 918 bp

**RefSeq ORF:** 597 bp

**Locus ID:** 84607

**UniProt ID:** [O88582](#)

**Cytogenetics:** 7q13

**MW:** 22.4 kDa

**Gene Summary:** suppressor of cytokine signalling [RGD, Feb 2006]