

## Product datasheet for **RR209826L4V**

### Zbtb25 (NM\_199496) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zbtb25 (NM_199496) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zbtb25
Synonyms:	KUP; Zfp50; Znf46
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199496
ORF Size:	1314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209826).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199496.1</a> , <a href="#">NP_955790.1</a>
RefSeq Size:	2633 bp
RefSeq ORF:	1317 bp
Locus ID:	314245
Cytogenetics:	6q24
Gene Summary:	putative transcription factor; human homolog is a DNA binding protein with two distantly spaced zinc fingers [RGD, Feb 2006]



[View online »](#)