

Product datasheet for **RR209833**

Rpl13a (NM_173340) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl13a (NM_173340) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpl13a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209833 representing NM_173340
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAGGGCAGGTTCTAGTATTGGATGGCCGGGCCATCTTCTGGCCGCCTGGCGCCATTGTGG
CCAAGCAGGTACTGCTGGGCCGAAAGTGGTGGTTGTACGCTGTGAGGGCATCAACATTTCTGAAATTT
CTACAGAAACAAGTTAAAGTATCTGGCCTTCTCCGAAAGCGGATGAACACCAACCCGCTCGAGGCCCC
TACCATTCCGAGCCCCAAGCCGATTTTTGGCGCACTGTGCGAGGCATGCTGCCGACAAAGACAAAA
GAGGCCAGGCTGCCCTGGAACGCCTCAAGGTGTTGGATGGATCCCTCCACCCTATGACAAGAAAAAGCG
GATGGTGGTCCCTGCTGCCCTCAAGGTTGTGCGGCTGAAGCCTACCAGAAAGTTGCTTACCTGGGGCGT
CTGGCTCATGAGGTGGGTGGAAGTACCAGGCAGTGACAGCTACTCTGGAGGAGAAACGGAAGGAAAAAGG
CAAAGATCCATTACCGAAGAAGAAGCAGCTCTTGAGGCTAAGGAAACAGGCAGAAAAGAATGTGGAGAA
GAAAATCTGCAAGTTCACAGAGGTCCTCAAGCCAATGGACTCTTGTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR209833 representing NM_173340
Red=Cloning site Green=Tags(s)

MAEQVLVLDGRGHLGRLAIVAKQVLLGRKVVVVRCEGINISGNFYRNKLYLAFLRKRMTNPSRGP
YHFRAPSRIFWRTVRGMLPHKTKRGQAALERLKVLDGIPPPYDKKKRMVPAALKVVRLKPTRKFAYLGR
LAHEVGWKYQAVTATLEEKRKEKAKIHYRKKQLLRKQAEKNVEKKICKFTEVLKTNGLLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_173340

ORF Size: 609 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173340.2</u> , <u>NP_775462.2</u>
RefSeq Size:	690 bp
RefSeq ORF:	612 bp
Locus ID:	317646
Cytogenetics:	1q22
MW:	23.4 kDa
Gene Summary:	60S ribosomal subunit protein; mouse homolog is mutated in the tum-transplantation antigen P198 [RGD, Feb 2006]