

Product datasheet for **RR209860L3V**

Plpp1 (NM_022538) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plpp1 (NM_022538) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Plpp1
Synonyms:	Ppap2a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022538
ORF Size:	846 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209860).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022538.2 , NP_071983.2
RefSeq Size:	1290 bp
RefSeq ORF:	849 bp
Locus ID:	64369
Cytogenetics:	2q14
Gene Summary:	hydrolyzes phosphatidic acid and lysophosphatidic acid; may mediate surfactant phospholipid secretion in alveolar type II cells, and cell growth during lung development and injury [RGD, Feb 2006]



[View online »](#)