

Product datasheet for **RR209872L3V**

Slc4a1 (NM_012651) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc4a1 (NM_012651) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc4a1
Synonyms:	MGC93554
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012651
ORF Size:	2784 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012651.2 , NP_036783.2
RefSeq Size:	4508 bp
RefSeq ORF:	2787 bp
Locus ID:	24779
UniProt ID:	P23562
Cytogenetics:	10q32.1
Gene Summary:	a transporter; functions as a chloride/carbonate exchanger and binds to components of the cytoskeleton [RGD, Feb 2006]



[View online »](#)