

Product datasheet for **RR209903**

Prmt2 (NM_001025144) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prmt2 (NM_001025144) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Prmt2
Synonyms:	Hrmt111
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR209903 representing NM_001025144
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCACCAGCTGAAGATCCAGCCATGAGTCACAGGTCACCCAGCTCCGGAGGAAGACCCTGTGG
 ACTATGGCTGTGAGATGCAGCTCCTTCAGGACGGGCACAGTTGCAGCTGCAGCTGCAGCCCAGGAGTT
 CGTAGCCATTGCAGATTACACCGCCACCGACGAGACACAGCTCAGCTTTCTGAGAGGAGAAAAAATTCTG
 ATCCTGAGGCAAACCACGGCGGACTGGTGGTGGGGCGAACGTGCGGGCTGCTGCGGGTACATCCCAGCAA
 ACCACCTTGGGAAGCAGGTGGAGGAGTACGACCCGGAGGACACCTGGCAAGATGAGGAGTACTTCGACAG
 CTATGGGACTCTGAACTTCACCTGGAATGTTGGCTGACCAGCCTCGGACAACCAAGTACCACAGTGTCT
 ATCCTGCAGAACAGGAGTCTCTGAAGGACAAGTCTACTGGATGTGGGCTGCGGCACCGGCATCATCA
 GCCTTTCTGCGCTCATCATGCCGGCCCAAGGAGTGTACGCGGTAGAAGCCAGTGACATGGCCCAGCA
 CACGGGCCAGCTGGTCTGCAGAATGGCTTTGCTGACACAATCACTGTGTTTCAGCAGAAGGTGGAGGAC
 GTGGTGTGGCCGAGAAGGTGGACGTGCTGGTGTGCGAGTGGATGGGGACCTGTCTCTTGTGTTGAGTTCA
 TGATCGAGTCCATCCTGTACGCCCGGATGCCTGGTTAAAGGAAGATGGCATCATTTGGCCAACCCAGGC
 GGCAGTGCAGTGGTACCCTGCAGCGCTGAGAAAGACTACCACAGCAAGGTGCTCTTCTGGGACAACCGG
 TACGAGTTCAACCTCAGTGCAGTCAAGTCGTTAGCAATCAAGGAGTCTTTTCAAGGCCAAAGTCTAATC
 ATATCTGAAGCCAGAAGACTGCCTCTCTGAGCCATGCACAATATTGCAGTTGGACATGAGAAGTGTGCA
 GGTCTCTGACCTGGAGACCATGCGTGGCGAGCTGCGATTTGACATCCAAAAGGCTGGGACCCTGCATGGC
 TTCACAGCTTGGTTCAGTGTACATTTCCAGAGCCTGGAGGAGGGCCAGCCTCAGCAGGTTCTCAGCACAG
 GACCGCTGCACCCACCACACACTGGAAGCAGACCTTGTATGATGGACGACCCAGTCCCAGTCCACAC
 AGGAGACGTGGTCACGGGTTGCGTGTACAGAGAAAATCCGGTGTGGAGAAGGCACATGTCCGTCTGT
 CTGAGCTGGTGTGCTACTTCCGAGTTGATCCCAGTCTCAAAGGGCTGGAGAGAAAGTCTTCCGATCT
 GGAGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR209903 representing NM_001025144
 Red=Cloning site Green=Tags(s)


MEAPAEDPSHESQVTPAPEEDPVDYGCEMQLLQDGAQLQLQPEEFVAIADYTATDETQLSFLRGEKIL
 ILRQTTADWWWGERAGCCGYIPANHLGKQVEEYDPEDTWQDEEYFDSYGTLLKHLLEMLADQPRTTKYHSV
 ILQNKESLKDKVILDVCGGTGIIISLFAHHPKAVYAVEASDMAQHTGQLVLQNGFADTITVFQKQVED
 VVLPEKVDVLVSEWMTCLLFEFMIESILYARDAWLKEDGIWPTTAAHLVPCSAEKDYHRSKVLFWDNA
 YEFNLSALKSLAIKEFFSRPKSNHILKPEDCLSEPCTILQLDMRTVQVSDLETMRGELRFDIQKAGTLHG
 FTAWFSVHFQSLLEEQPQQVLSGTLPHPTTHWKQTLFMMDDPVPVHTGDVVTGSVVLQRNPVWRRHMSVC
 LSWVTSAVDPTSQRAGEKVFPIWR

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001025144

ORF Size: 1335 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025144.1](#), [NP_001020315.1](#)

RefSeq Size: 2014 bp

RefSeq ORF: 1338 bp

Locus ID: 499420

Cytogenetics: 20p12

MW: 50.5 kDa