

Product datasheet for **RR209939L3V**

Lgals3bp (NM_139096) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lgals3bp (NM_139096) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lgals3bp
Synonyms:	CyCAP; MAC2BP; mama; Ppicap
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_139096
ORF Size:	1722 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139096.1 , NP_620796.1
RefSeq Size:	2151 bp
RefSeq ORF:	1725 bp
Locus ID:	245955
UniProt ID:	O70513
Cytogenetics:	10q32.3
Gene Summary:	displays differential expression in a fibroblast cell line transformed by human T-cell leukemia virus type 1 Tax protein [RGD, Feb 2006]



[View online »](#)