

Product datasheet for **RR209963L3V**

Nqo2 (NM_001004214) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nqo2 (NM_001004214) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nqo2
Synonyms:	MGC94180
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004214
ORF Size:	693 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209963).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004214.1 , NP_001004214.1
RefSeq Size:	1235 bp
RefSeq ORF:	696 bp
Locus ID:	291084
UniProt ID:	Q6AY80
Cytogenetics:	17p12
Gene Summary:	mouse homolog catalyzes the reduction of quinones and quinone derivatives; may protect against development of myelogenous hyperplasia , may contribute to menadione mediated hepatic toxicity [RGD, Feb 2006]



[View online »](#)