

Product datasheet for RR209966L4

Hnrnpu (NM_057139) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hnrnpu (NM_057139) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Hnrnpu
Synonyms:	hnRNP U; Hnrpu; SAF-A; SN1; SP120
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209966).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_057139
ORF Size:	2394 bp



[View online »](#)

