

Product datasheet for **RR209975L4V**

Etfa (NM_001009668) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Etfa (NM_001009668) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Etfa
Synonyms:	ETF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009668
ORF Size:	999 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209975).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009668.1 , NP_001009668.1
RefSeq Size:	1426 bp
RefSeq ORF:	1002 bp
Locus ID:	300726
UniProt ID:	P13803
Cytogenetics:	8q24
Gene Summary:	specific electron acceptor for several dehydrogenases, including acyl-coA dehydrogenases, glutaryl-coA and sarcosine dehydrogenase; transfers electrons to the main mitochondrial respiratory chain [RGD, Feb 2006]



[View online »](#)