

Product datasheet for **RR209978**

Ccm2 (NM_001126275) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccm2 (NM_001126275) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ccm2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR209978 representing NM_001126275
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGAGGAGGGCAAGAAGGGCAAGAAGCCTGGAATTGTGTCGCCGTTAAACGAGTATTCCTAAAAG
GTGAAAAGAGTAGAGATAAGAAAGCCCATGAGAAGGTGACAGAGAGGCCCCCTCTACACACTGTGGTGCT
GGCCTTGCTGAGCGAGTCGAGCCAGACAGACTGCTGAGCGACTATATTGAGAAGGAGGTGAAGTATTTA
GGCAGTTAACATCCATTCTGGCTACCTGAACCCTCCAGTAGGACGGAAATCCTGCATTCATAGACA
AGGCAAAGCGGTCCACCAACTCCCTGGCACCTGACTCAGGAGCATGATGCTGTGCTCAGTCTGTCCGC
CTACAATGTCAAGCTGGCTGGAGGGATGGGGAGGACATTATCCTCAGGGTGCCCATCCACGACATCGCT
GCTGTCTCCTATGTCGAGATGATGCTGCACACCTGGTGGTCTGAAGACAGCCAGGACCCAGGCATCT
CTCCAGCCAAAGTCTGTGTGCGGAAGTTCTAGAGGCCCTCAGCGCAGGCTCCTTGTGAGAAAGTGCAGT
GGGGCCAGTAGAGGCTGTTGCCTGGTTCATCATGGCCACAGAGCAAGGTGCGCCGCTGAAGAACTGTGC
TCCCTGCTCAGCCAGTCTTCCAGATCGTCTACACAGAGTCCACCATCGACTTTCTGGACCGAGCAATAT
TCGACGGGGCCTCCACCCTACCCACCCTGTCAGTGCACAGCGATGACTCTTCCACGAAAGTGGACAT
GAAGGACAGTTATGACGCCGACGCCAGCACCTTCTGCTTCCCGGACTCTGGGGATGTGGGGGCTTGCCG
CCCTTGCCCTTCTGCATGCAGGCATCCCCCACAGCAAGACGGTACGCGAGAGTGAGCTGAGTACCAGCG
CCACGGAGCTGCTGCAGGACTACATGCTCACGTTACGTACCAAGTGTGTCATCACAGGAGATCCAGCAGTT
TGCAGCCCTGCTGCATGAATACCGCAACGGGGCCTCCATCCATGAGTTCTGCATCAGCCTGCGGCAGCTC
TATGGAGACAGCCGAAGTTCTGCTGCTTGGTCTGAGGCCCTTATACCTGAGAAGGACAGTCACTGACT
TTGAAAATTCTGGAGACCATTGGCGTAAAGACGGCCGTGGCATCATCACTGACAGCTTTGGCAGGCA
TCGCCGTGCTCTGAGCACCACTCCACATCCACCACCAATGGGAACAGGACCACAGGCAGCCCTGATGAC
CGCTCAGCGCCCTCAGAGGGGGATGAGTGGGACCGCATGATCTCAGACATCAGCAGTGATATTGAAGCAC
TAGGCTGCAGCATGGACCAGGACTCAGCG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA




Protein Sequence: >RR209978 representing NM_001126275
 Red=Cloning site Green=Tags(s)

MEEEGKKGKPGIVSPFKRVFLKGEKSRDKKAHEKVTERPLHTVVLALPERVEPDRLLSDYIEKEVKYL
 GQLTSSIPGYLNPSRTEILHFIDKAKRSHQLPGHLTQEHDVLSL SAYNVKLAWRDGEDILRVPIHDIA
 AVSYVRDAAHLVVLKTAQDPGISPSQSLCAEGSRGLSAGLSSEAVGPVEACCLVIMATESKVAEEELC
 SLLSQVFQIVYTESTIDFLDRAIFDGASTPTHHLSVHSDSSTKVKMDKSYDADASTFCFPDSGDVGGGLP
 PLPFCMQASPHTVSESELSTSATELLQDYMLLRTKLSSQEIQQAALLHEYRNGASIHFCISLRQL
 YGDSRKFLLLGLRPFIPKDSQHFENFLETIGVKDGRGIIITDSFGRHRRALSTTSTTTNGNRTTGSPDD
 RSAPSEGDEWDRMISDISSDIEALGCSMDQDSA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001126275
ORF Size:	1359 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001126275.1 , NP_001119747.1
RefSeq Size:	1807 bp
RefSeq ORF:	1362 bp
Locus ID:	305505
Cytogenetics:	14q21
MW:	49.8 kDa