

## Product datasheet for **RR209980**

### **Ebf2 (NM\_001108383) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Ebf2 (NM_001108383) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ebf2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RR209980 representing NM\_001108383  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTTTGGGATTCAAGATACGCTAGGAAGAGACCAGTCTGAAAGAAAAGTCGCTGGTGCCGAGATGG  
 ATTCGGTCAGGTCTGGGTCCGGAATGTCGGCGTGGTGGATGCCAATGTCGCGGCTCAGAGTGCCGTTGC  
 CCTGTCTCGGGCCACTTTGAAAAGCAACCTCCTTCTAACTTGAGGAAATCCAACCTCTTCCATTTTGT  
 CTGGCGCTCTACGACAGGCAGGGACAGCCGGTGGAGATTGAGCGCACAGCCTTCGTGGACTTTGTGGAGA  
 ATGACAAAAGCAAGGCAACGAGAAAACCAACAATGGCACTACTACAAGTTACAGCTCCTTTACAGCAA  
 TGGTGTCCGCACAGAGCAGGACCTCTATGTCAGGCTCATAGACTCAGTCACAAAACAGCCTATTGCCTAT  
 GAGGGACAGAATAAGAATCCAGAAATGTCAGAGTTCTGCTGACACAGAGGTGATGTGCAGCCGATGCT  
 GCGAGAAAAGAGCTGTGGGAACCGGAATGAGACCCCGTCAGACCCCGTCATCATCGACAGATTCTTTTT  
 AAAATTTTCTGAAGTGAATCAGAATTGTTTAAAACAGCAGGAAACCAAGGGACATGAGACGGTTT  
 CAGGTCGTGTTGTCTACAACGGTCAATGTGGATGGACATGTGCTGGCTGTTTCTGACAACATGTTTGTTC  
 ATAACAACCTAAGCATGGAAGGAGAGCGAGAAGGCTTGACCCATCAGAAGCTACACCCTGCATCAAAGC  
 CATCAGTCCAAGTGAAGGCTGGACCACGGGAGGGGCCATGGTCATCATCATTGGTGACAACCTTTTTGTAT  
 GGTCTCCAAGTTGATTTGGTACCATGCTCGTATGGAGTGAGCTTATAACCCCTCATGCCATCAGAGTAC  
 AGACCCCTCCCCGGCATATCCCTGGAGTTGTGGAAGTGACACTATCTTATAAATCCAAACAGTTCTGCAA  
 AGGAGCCCCAGGCAGGTTTCATCTACACAGCCTTAAATGAACCCACCATAGACTATGGCTTCCAAAGACTG  
 CAGAAAGTCATCCAAGGCACCCAGGAGACCTGAGAGGCTAGCTAAGGAAATGCTATTGAAAAGAGCTG  
 CAGACCTGGTGAAGCCCTCTACGGCACACCCCAACAATAACAGGACATCATTTTAAAGCGTGTGCAGA  
 CATTGCAGAAGCCCTGTACAGTGTCCCGGAACCCAGCCAGATTCCAGCTCTCTCCAGCTCCCCAGCT  
 CATAGCGCATGATGGGTATCAACTCTTACGGCAGCCAGCTGGGGTTCAGCATTTCAGAGTCTACACAAG  
 GAAATAACCAAGGATACATCCGCAACACAAGCAGCATCTCCCTCGGGGATACTCTCCAGCTCCACGCC  
 CCAGCAGTCTAATTACAGCACCTCCAGTAACAGCATGAATGGCTACAGCAACGTTCTATGGCCAACCTG  
 GGTGTGCCAGGATCACCGGTTTTCTGAACGGCTCACCTACAGGCTCCCCATATGGAATCATGTCTTCAA  
 GCCCCACCGTGGGATCTCCAGCACATCTCCATCTCCCGTTCTCTCTCAGTTTTCTCTGTGTCAA  
 ACAGAAGAGTGCCTTTGCCCTGTATCAGGCCCAAGGCTCCCCTTCGCTGCCTGCTCCAGTGGCAAT  
 GGAATGGATTACAGCCATGACGGGACTTGTGTACCCCGATG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RR209980 representing NM\_001108383  
 Red=Cloning site Green=Tags(s)


MFGIQDTLGRGPVLKEKSLGAEMDSVRSWVRNVGVVDANVAAQSGVALSRAHFQKPPSNLRSNFFHFV  
 LALYDRQGPVEIERTAFVDFVENDKEQNEKTNNGTHYKQLLYSNGVRTEQDLVYRLIDSVTQKPIAY  
 EGQNKNPMECRVLLTHEVMCSRCEKKSNGNRNETPSDPVIIDRFLLKFLKCNQNLKTAGNPRDMRRF  
 QVVLSTTVNVDGHVLAUSDNMFVHNSKHGRRARRLDPSEATPCIKAIISPSEGWTGGAMVIIIGDNFFD  
 GLQVVFVGTMLVWSELITPHAIRVQTPPRHIPGVVEVTLVSYKSKQFCKGAPGRFIYALNEPTIDYGFQRL  
 QKVIPRHPGDPERLAKEMLLKRAADLVEALYGTPHNNQDIILKRAADIAEALYSVPRNPSQIPALSSSPA  
 HSGMMGINSYGSQLVSIESTQGNQGYIRNTSSISPRGYSSSTPQQSNYSTSSNSMNGYSNVPMANL  
 GVPGSPGLNGSPTGSPYIMSSSPTVGSSTSSILPFSSVFPVAVKQKSAFAPVIRPQGSPPACSSGN  
 GNGFRAMTGLVPPPM

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM\_001108383

<b>ORF Size:</b>	1725 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001108383.1</a> , <a href="#">NP_001101853.1</a>
<b>RefSeq Size:</b>	5439 bp
<b>RefSeq ORF:</b>	1728 bp
<b>Locus ID:</b>	361060
<b>Cytogenetics:</b>	15p12
<b>MW:</b>	62.6 kDa