

## Product datasheet for **RR209984L3V**

### Slc25a47 (NM\_001001509) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a47 (NM_001001509) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a47
Synonyms:	Hdmcp; Slc24a47
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001001509
ORF Size:	930 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209984).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001001509.1</a> , <a href="#">NP_001001509.1</a>
RefSeq Size:	1550 bp
RefSeq ORF:	933 bp
Locus ID:	299316
UniProt ID:	<a href="#">Q6J329</a>
Cytogenetics:	6q32
Gene Summary:	Uncoupling protein which may catalyze the physiological 'proton leak' in liver. [UniProtKB/Swiss-Prot Function]



[View online »](#)