

Product datasheet for **RR209999L4V**

Slc10a2 (NM_017222) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc10a2 (NM_017222) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc10a2
Synonyms:	ASBT; IBAT; ISBAT; ISBT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017222
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR209999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017222.2 , NP_058918.1
RefSeq Size:	4269 bp
RefSeq ORF:	1047 bp
Locus ID:	29500
UniProt ID:	Q62633
Cytogenetics:	16q12.5
Gene Summary:	mediates sodium-dependent transport of bile acids; plays a role in bile acid circulation and reabsorption; facilitates cholesterol catabolism [RGD, Feb 2006]



[View online »](#)