

Product datasheet for **RR210074L3V**

Gtf2h4 (NM_212501) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtf2h4 (NM_212501) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gtf2h4
Synonyms:	MGC156549
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_212501
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210074).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_212501.2 , NP_997666.1
RefSeq Size:	1692 bp
RefSeq ORF:	1392 bp
Locus ID:	294236
Cytogenetics:	20p12
Gene Summary:	human homolog is a subunit of the transcription factor TFIIF complex and plays roles in transcription and DNA repair [RGD, Feb 2006]



[View online »](#)