

Product datasheet for **RR210078L4V**

Trdmt1 (NM_001031643) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trdmt1 (NM_001031643) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trdmt1
Synonyms:	Dnmt2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001031643
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210078).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001031643.1 , NP_001026813.1
RefSeq Size:	1536 bp
RefSeq ORF:	1176 bp
Locus ID:	291324
UniProt ID:	Q4G073
Cytogenetics:	17q12.3
Gene Summary:	Specifically methylates cytosine 38 in the anticodon loop of tRNA(Asp).[UniProtKB/Swiss-Prot Function]



[View online »](#)