

Product datasheet for **RR210088L3V**

Birc2 (NM_021752) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: Birc2

Synonyms: Api2; rIAP1

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_021752

ORF Size: 1767 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RR210088).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_021752.2](#), [NP_068520.2](#)

RefSeq Size: 2907 bp

RefSeq ORF: 1770 bp

Locus ID: 60371

Cytogenetics: 8q11

Gene Summary: may play a role in regulation of apoptosis [RGD, Feb 2006]

