

Product datasheet for **RR210097**

Pot1 (NM_001024322) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pot1 (NM_001024322) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pot1
Synonyms:	Pot1a
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR210097 representing NM_001024322
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCTTTGGTTCGGTAACTCCCTATACATATACACCCCTGAATCTACTTAAGGAAGTGCAGTTGCCA
 ATGTCTATGGTGTGTGAAATCTTTAAGCCTCCATATCTCAGCAAAGGAAGTATTATTGCTCAGTTGT
 AACAAATTGTGGACCCAACAAATGTGAAGTTAACGTGCATGCTCTTTAGTGAAACTATGAAGCCCTTCCG
 ATCATTACAAAGTTGGAGACATTGTTTCGCTTCCATAGGCTGAAGATCCAAGTATATAAAAAATGAGCTTC
 AGGGAATCAACAGCTCTGGGTTTCGCTCTTTGACATTTGAAGGAACATTAGGGACACCTGTCACAGCAGC
 TACTTCAAGCAAAGTTTTAGCTTCACTTCTCAGGACCAAAAAATGGTAGAAGCTTTGCGAGTTTGGGCA
 TCCAAGCACATTGCAGCTTCTTCAAGTTTAGTGCAGCTGTGTGACGCTCAGCCCCGCAATATTATGACC
 TGACCTGTGAGCTTCTGGGTAAGCAGAAGTGGATGGAACAGCCTTTCTTCTAAAGGTATGGGACGGAAC
 CAGGACTATGTTCCCTCTTGGAGAGTCCATACAAAATCTTACCTTTGAAGGTGACTTGAGTCACATT
 CAGCGGCTACAGAGTCTGGTGATAGACATTCTAGTCTATGACAATCATGTTCAAGTGGCAAAATCCCTCA
 AGATTGGAAGCTTTATTGAATCTATAGCTCCATACAAAATCTCAGTCTGTTGATTTCAGAAACCAGCAA
 GTCATTGTTAAGGCTGGAGTTTCATCTCCATGGAGGGACAGTTATGGTGGGGGATTAGGGTTTTGCCA
 GATGAAAGCTCGTGTGTGGATCAGCTGAAAAAGATTTAGAATGTGCAGATTCGGCAGTCACTCAACCT
 CAAATGGTATCTGTCAATCAAAAAATGAGGATAGCTTTGTACTTTCAAGCTCTGCATCAGGCTCAGTATC
 CCTGTATGAGGAAGAACGATGTCAGCAAGTCTCTTCTACAATACTCACAGATCACCAGCACTTAGAGAAA
 ACTCCATTGTGTGCCGTTTTGAAGCAAAGGCTCCGACGCAATATCGGGTCCGTGCCAAGCTGAGGTCGT
 ATTTGCCAGGAGACTTCCCAGTCTGTTAAACTGATCTGTCCGAGTCCATTACTACTCAAGAAGTTCC
 ACATGAGGACAGTTTGGATAAACTTTTGAAGATGCTGCAACTAAAGCGCCAGACATCAGGCTGCAAGGT
 ACATCTTTGTATGAATCAAAGATCTGGACCACCGAAGGCCAGGGAGGACGGCAGGTGCTGTTCACTTTG
 TGAAGAAGGATGGTGTCTCCACTTTCAAGTGAATGTCTAATTTTATAGAAGGAGGAAGACTCTGTGA
 AGTGTCCAAGCTGTCCAGCAAGTCCACAGTGTGATCCCTGTGAGGTCCGGCCCTGAGGATCTGGAGCTC
 CTCGACATTTAGCCCCATTCTCATACAGGAAAGCTGCATCACTATGGATGTAACAGTGTTCAGCT
 TGAACCTATACAAAACCTAAATTCCTTCTTAATAAAGGACCATGGATTCCTTCTCCGTGGCAAAAGA
 GATGCTGATAGTATCCTTTGGCAACCCAAACCGACAACCAGATCTTTTTGTTAGACTGGAATTTTATA
 AATTCCTGGGCATTGTGCTCTGCAGTATGTGTTGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR210097 representing NM_001024322
 Red=Cloning site Green=Tags(s)


MSLVPTPYTYTPLNLLKEGAVANVYGVVKKFKPPYLSKGTDYCSVVTIVDPTNVKLTCLFSGNYEALP
 ILYKVDIVRFHRLKIQVYKNEQLQINSSGFASLTFEGLTGTPVTARTSSKVFSTSDQDKMVEALRWVA
 SKHIAASSSLVQLCDAQPLQYYDLTCQLLGAEVDGTAFLKLVWDGTRTMFPSWRVSIQNLTFEGDLSHI
 QRLQSLVIDILVYDNHVQVAKSLKIGSFIRIYSLHTKLQSVSETSKSLLRLEFHLHGGTSYGRGIRVLP
 DESSVDQLKKDLECADSAVTQPSNGICQSKNEDSFVLSASGSVSLYEEERCQVSSITLDHQHLEK
 TPLCAVLKQKAPQYRVRALRSYLPRLPQSVKLCPCQCHSLQEVPHEDSLDKLLQDAATKAPDIRLQG
 TSLYESKIWTTEGQGRQVAVHFVKKDGVLPSSSECLILEGGRLCEVSKLSSKFHVSIPVRSVPEDLEL
 LDISAPFLIQGLHHYGCKQCSSLKPIQNLNSLLNKGWPWIPSSVAKEMLIVSFGQPKPTTRFFLLDWNFI
 NSWALCLCSMCLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001024322
 ORF Size: 1719 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001024322.1 , NP_001019493.1
RefSeq Size:	2730 bp
RefSeq ORF:	1722 bp
Locus ID:	500054
Cytogenetics:	4q22
MW:	63.8 kDa