

## Product datasheet for **RR210099**

### Fam3b (NM\_001107102) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fam3b (NM_001107102) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fam3b
Synonyms:	RGD1311662
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR210099 representing NM_001107102 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGTCCAGTTACAGCAGGTATCTTCAAGGCACTACTGTTTATTTCTCCTCCCTGTGCGCCTGGTATT  
CTGGGTACCTGCTCGCGGAGCTCATTCTGACGTGCCCTGTCCAACACTCTCTACAACATCCGAAGCAT  
TGGGGAGAGACCTGTTCTCAAAGCCCCAGTCCCCAAAAGACAAAATGTGACCATTGGTCCCATGTCCT  
CCTGACACCTATGCCTACCGGCTGCTCAGTGGTGGCCTAGACAAGTATGCCAAGATCTGCTTTGAGG  
ATGAAATGCTATTAGGAGAGAAGACGGGAATGTGGCAAGAGGGATAAACATCGCTGTTGTTGACTATGA  
GACGGGAAAAGTAATATCAACAAAGTACTTTGACATGTACGAAGGTGATAACTCTGAGCCAATGACTAAG  
TTTATTCAGAGCATCCCTTCGAAATCCCTGCTGTTTCATGGTGACTCATGATGACGGAAGCCACAACTGC  
AGGCTCAAGCAAAGGATGCCATAGAAGCCCTTGAAGCAAAGAAATCAAGAACATAAAGTTCAGATCAAG  
CTGGGTGTTTGTGTCAGCAAAGGCTTTGTGCTCCCTTCAGAAAGGAGAGAGAAAAATCAACCCTCAG  
ATCAATCCAGGAACAGATACTCCGGCTGCCAGCAGAGATCCAGATCGAAGGATGCATACCCAAAGGGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR210099 representing NM\_001107102  
 Red=Cloning site Green=Tags(s)

MRPVTAGIFKALLFIFSSLCAWYSGYLLAELIPDVPLSNTLYNIRSIGERPVLKAPVPKRQKCDHWSPCP  
 PDTYAYRLLSGGGLDKYAKICFEDEMLLGEKTNVARGINIAVVDYETGKVIKSTKYFDMYEGDNSEPMTK  
 FIQSIPSKSLLFMVTHDDGSHKLQAQAKDAIEALGSKEIKNIKFRSSWVFVAAKGFVLPSSERREKKSTTQ  
 INPGTDTPAGQQRSRKDAYPKG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001107102

**ORF Size:** 699 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001107102.1</a> , <a href="#">NP_001100572.1</a>
<b>RefSeq Size:</b>	952 bp
<b>RefSeq ORF:</b>	702 bp
<b>Locus ID:</b>	304050
<b>Cytogenetics:</b>	11q12
<b>MW:</b>	25.9 kDa