

Product datasheet for **RR210106L3V**

Mustn1 (NM_181368) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mustn1 (NM_181368) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Mustn1
Synonyms:	Mustang
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181368
ORF Size:	246 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210106).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181368.3 , NP_852033.1
RefSeq Size:	1097 bp
RefSeq ORF:	249 bp
Locus ID:	290553
UniProt ID:	Q80XX4
Cytogenetics:	16p16
Gene Summary:	a nuclear protein thought to be involved in bone regeneration after injury and bone formation during embryogenesis [RGD, Feb 2006]


[View online »](#)