

Product datasheet for **RR210116L4V**

Lrrtm3 (NM_001106387) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrtm3 (NM_001106387) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrtm3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001106387
ORF Size:	1746 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210116).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001106387.1 , NP_001099857.1
RefSeq Size:	3837 bp
RefSeq ORF:	1749 bp
Locus ID:	294380
UniProt ID:	D3ZAL8
Cytogenetics:	20p11
Gene Summary:	May play a role in the development and maintenance of the nervous system (By similarity). Exhibits a limited synaptogenic activity in vitro, restricted to excitatory presynaptic differentiation.[UniProtKB/Swiss-Prot Function]


[View online »](#)