

Product datasheet for **RR210140**

LOC501038 (NM_001047965) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: LOC501038 (NM_001047965) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: LOC501038
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR210140 representing NM_001047965
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTAACGTCAGTACTAACAATGAACTGCTACAAGGTTTATAGCGTCCACACCATTAAGAGTGGTTCAT
 TCAAAGTGCCTAAACAGATAGTGTGTAAAGTCAAAGCAGGACCACACATACACTCTGCATGGAGAA
 GCTGCTCCTAGGCTTCCCTGGACTGTGTAAGATCCTGCCAGTGGAGTAACTGCCATGCCTTTGTGAAG
 TCAGCTGTCTCAAAGTACGAACAGAAACCACTGACAGCACTACATTCCAGGTCTTGATCCTCTACAACA
 TAGAAAAATGTGTAGGAACTTCCAGAAGGGTGGACATCGTGTGGATGTTATCAGTATAAAACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210140 representing NM_001047965
 Red=Cloning site Green=Tags(s)

MVTSVLTMNCYKVYSVHTIKSGSFKVPKTDVVKSKQDHTHTLCMEKLLLGFPLCKDPASGVNCHAFVK
 SAVSKYEQKPLTALHFQVLILYNIIEKCVGNFQKGGHRVDVYQYKT

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM 001047965

ORF Size: 345 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM 001047965.1, NP 001041430.1

RefSeq Size: 394 bp

RefSeq ORF:	348 bp
Locus ID:	501038
Cytogenetics:	8q31
MW:	12.9 kDa
Gene Summary:	DISCONTINUED: This record has been withdrawn by NCBI Staff. This record was defined by AY325184.1, which is not sufficient evidence to define a distinct gene.

Product images:



Circular map for RR210140