

Product datasheet for **RR210140L4V**

LOC501038 (NM_001047965) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LOC501038 (NM_001047965) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	LOC501038
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001047965
ORF Size:	345 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210140).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001047965.1 , NP_001041430.1
RefSeq Size:	394 bp
RefSeq ORF:	348 bp
Locus ID:	501038
Cytogenetics:	8q31
Gene Summary:	DISCONTINUED: This record has been withdrawn by NCBI Staff. This record was defined by AY325184.1, which is not sufficient evidence to define a distinct gene.


[View online »](#)