

Product datasheet for **RR210158L4V**

Serpina1 (NM_022519) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serpina1 (NM_022519) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Serpina1
Synonyms:	AAT; Pi; Spi1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022519
ORF Size:	1233 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210158).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022519.2 , NP_071964.2
RefSeq Size:	1603 bp
RefSeq ORF:	1236 bp
Locus ID:	24648
UniProt ID:	P17475
Cytogenetics:	6q32
Gene Summary:	Inhibitor of serine proteases. The primary target is elastase, but also has a moderate affinity for plasmin and thrombin.[UniProtKB/Swiss-Prot Function]



[View online »](#)