

Product datasheet for RR210227

Polr1e (NM_001107938) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Polr1e (NM_001107938) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Polr1e
Synonyms: Praf1; RGD1565773
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210227 representing NM_001107938
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTACGCCGATTCTCCCGGATGGATGGCCAGACTGGAGACACGGAGACAGAGGTGCTACAGAGTG
CAAGATGGCAGTACTGTGGAGAGCCCGACGACAGGCAGAAAGCTGTGCTGGTCCAATTCTCCAATGGAA
GCTGCAGAATCCAGGCAACATGCGCTTACCTTGTACAACAGCACCGACTTGGTAAACCCAGGCAGAAG
TGCCATCGGATCCTGGCCGCTGAGACAGACAGGCTCTCCTATGTTGGAACAATTTGGGACTGGCGCC
TCAAAATGCAACACTCTTTGCAGACACTTTGTGGCATTCTGAATAAGACCTCTGGCCAAATGGAAGTCTA
TGATGCGGAAGTATTCAACATGCAGCCCCTGTTTGTGAAGACTCTATTGAGCATGGGCCTCCTCTGGAG
AATCAGAAACAAGACCTTCAGGGACAAGTTGGATTCTTGCATTGAAGCCTTTGGAAGCACCAAGCAGAAGC
GGTCCCTGAACCTCCCGGAGAATGAACAAAGTTGGCAGCGAGTCTTTGAATTTACAGTAGCTAAAGCCG
AGAGAGTATCATTGACACGAAGGGCGTGAGTGTCTTGGTCAGCGATGCCATGCAAGATGACTTGCAAAAT
GACTCCCTACCTTCTCCCTGCCATGCTGATGCCACCAAGCCGAAGACGTGTATAGATTTGAAGACA
TCCTGTCCCGGCGAGTATGATGCTCTTGAAGCCCATCCGAAGCCTTCAGAAAGTCTATGCTGAAGA
CATACTGAAGATGGTGAAGAGAACAGCCACTGCTCCTTCATCATAGAAATGCTGAAGTCTTTGCCAGCA
GACGAGGTGCAGCGAAACCGCAAGCCCGAGCATTGTTTCTAGATGCCCTCCTCAGATTCAGGCTC
AGAAAAGTATAAAAGGGAAAAGTGCCTTAGGACCTGGGATCCCGCACATCATCAACACCAAGCTGCTGAA
GCAGTTACCTGCTTGACCTACAACAACGGCAGTTTACGGAACCTAATCTCGAGTTCTATGAAGGCCAAG
ATCACTGCATATGCAATCATCTTGTCTGCACATAAATAACTTCCAGATTGATCTGACGGTGTACAGA
GGGACTGAAGCTGAGTGAGAAAAGGATGATTGAGATAGCGAGAGCCATGAGGCTGAAGATCTCCAAAAG
AAAGGTGTCCTGGCAGATGGCAGGGAAGAGGATCATAGACTGGGCACTCTGTCCGTCCTGCTCCG
GCCAGACCTCAGACCGACAGTCAAGCGGAAAAAATGAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR210227 representing NM_001107938
 Red=Cloning site Green=Tags(s)


MATPDSPGMDGQTGDTETEVLQSARWQYCGEPDDRQKAVLVQFSNGKLQNPNGMRFTLYNSTDLVNPQRK
 CHRILAAETDRLSYVGNFNGTGALKCNTLCRHFVGI LNKTSQMEVYDAEVFNMQPLFAEDSIEHGPPLE
 NQNKTRFDKLDSCIEAFGSTKQKRSLNSRRMNKVGSESLNFTVAKAAESI IDTKGVSA LVSDAMQDDLQN
 DSLYLPPCHADATKPEDVYRFEDILSPA EYDALESPSEAFRKMVSEDI LKMVEENSHCSF I IEMLKSLPA
 DEVQRNRQARSIWFLDALLRFRAQKVIKGSALGPGIPI INTKLLKQFTCLTYNNGSLRNLISSSMKAK
 ITAYAILALHINNFDLTVLQRDLKLSEKRMIEIARAMRLKISKRKVSLADGREEDHRLGTL SVPLPP
 AQTSDRQSKRKKMS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107938

ORF Size:	1302 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107938.1 , NP_001101408.1
RefSeq Size:	1588 bp
RefSeq ORF:	1305 bp
Locus ID:	313245
Cytogenetics:	5q22
MW:	48.8 kDa