

## Product datasheet for **RR210237L4V**

### IL22ra2 (NM\_001003404) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL22ra2 (NM_001003404) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	IL22ra2
Synonyms:	Crf2-s1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001003404
ORF Size:	687 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210237).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001003404.1</a> , <a href="#">NP_001003404.1</a>
RefSeq Size:	2506 bp
RefSeq ORF:	690 bp
Locus ID:	444986
UniProt ID:	<a href="#">Q7TNI4</a>
Cytogenetics:	1p12
Gene Summary:	Receptor for IL22. Binds to IL22, prevents interaction with the functional IL-22R complex and blocks the activity of IL22 (in vitro). May play an important role as an IL22 antagonist in the regulation of inflammatory responses.[UniProtKB/Swiss-Prot Function]



[View online »](#)