

## Product datasheet for **RR210240L3V**

### Hes5 (NM\_024383) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hes5 (NM_024383) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hes5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024383
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210240).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024383.1</a> , <a href="#">NP_077359.1</a>
RefSeq Size:	592 bp
RefSeq ORF:	501 bp
Locus ID:	79225
UniProt ID:	<a href="#">Q03062</a>
Cytogenetics:	5q36
Gene Summary:	basic helix-loop-helix (HLH) protein; may function as negative regulator in neural development [RGD, Feb 2006]



[View online »](#)