

Product datasheet for RR210245L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Ago2 (NM_021597) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Ago2 (NM_021597) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Ago2 Synonyms: Eif2c2

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_021597 **ORF Size:** 2589 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RR210245).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 021597.1, NP 067608.1

 RefSeq Size:
 2626 bp

 RefSeq ORF:
 2592 bp

 Locus ID:
 59117

 UniProt ID:
 Q9QZ81

 Cytogenetics:
 7q34

Gene Summary: may have important roles in development and exert the effects on cell differentiation from

the level of intracellular membranes [RGD, Feb 2006]

