

Product datasheet for **RR210256L3V**

Kpnb1 (NM_017063) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kpnb1 (NM_017063) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kpnb1
Synonyms:	Impnb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017063
ORF Size:	2625 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210256).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017063.1 , NP_058759.1
RefSeq Size:	2991 bp
RefSeq ORF:	2628 bp
Locus ID:	24917
UniProt ID:	P52296
Cytogenetics:	10q31
Gene Summary:	subunit of complex necessary for nuclear import [RGD, Feb 2006]



[View online »](#)