

## Product datasheet for **RR210277L4V**

### Canx (NM\_172008) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Canx (NM_172008) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Canx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172008
ORF Size:	1773 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210277).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172008.2</a> , <a href="#">NP_742005.1</a>
RefSeq Size:	4986 bp
RefSeq ORF:	1776 bp
Locus ID:	29144
UniProt ID:	<a href="#">P35565</a>
Cytogenetics:	10q22
Gene Summary:	calcium-binding protein that may serve to chaperone proteins from the ER to outer cellular membrane [RGD, Feb 2006]



[View online »](#)